AiSEMI simple Temperature Measurement System Manual

With array-typed thermopile module sensor and integrated prime lens, the AiSEMI simple can measure surface skin temperature accurately when used correctly. Please install the AiSEMI simple on the optional Tripod or placed on the table. Power on the system with a 5V_1A USB charger (Alkaline Batteries are only for backup). Position the person being evaluated at 50cm distance directly facing the AiSEMI simple for 1-3 seconds when used. **Distance will affect the result of temperature measurement.**

Please keep the system in a room with no draft (movement of air). For best results, avoid direct sunlight and away from radiant heat or AC vent. Please make sure that there is no 940nm light source in front of this product, such as sunlight, TV screen light, or camera with night vision function, etc.

If the body temperature exceeds 37.5°C or 99.5°F, AiSEMI simple will turn on a red light with a warning sound.

***To obtain accurate temperature measurement, please turn on the AiSEMI simple 15 minutes before being used to warm it up so the temperatures between system and the room temperature balance out.**

*** Please place the AiSEMI simple at 130cm height. Avoid strong lighting and away from radiant heat or AC vent.**

*** Position a floor sign at a fixed distance 50cm from the AiSEMI simple system, directly facing the camera before measurement.**

Out-of-the-box inspection:

(A) AiSEMI simple Main Body x1

(B) Type-C to USB Cable x1

(C) Positioning Sticker x1

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(E) 5V_1A USB Adapter (Excluded)


Turn (A) over, lift the back cover, connect (B) Type-C power cable or alternatively install batteries (4* AA Alkaline Batteries) for use.

(A) can be placed on the table (at a height of about 110cm) or installed on the tripod (D) (at a height of about 130cm.)

Button and Function Description :

After powering on, please press **+** upward adjustment key to turn on the machine, the screen will appear, and the red LED will light up and turn off automatically. Press **+** upward adjustment key and hold for 6 seconds to shut down the machine.

After pressing **SET** mode setting key, these settings can be made in sequence:

Temperature Unit, Warning Temperature, High Temperature Warning, Multi-lingual Voice, Volume.

1. Temperature Unit : One-click **SET** to C/F Temperature Unit icon. Press **+** or **-** to select the unit of °C (Celsius) or °F (Fahrenheit). After setting, press **SET** again or wait 5 seconds for it to finish automatically.

2. Warning Temperature : Double-click **SET** to Warning Temperature icon. This is the warning value of high temperature. Default is 37.5°C (99.5°F). Press **+** or **-** to select the warning value. When the temperature exceeds this value, AiSEMI simple will turn on a red light and give off a warning voice. After setting, press **SET** again or wait 5 seconds for it to finish automatically.

3. High Temperature Warning : Three-click **SET** to High Temperature Warning icon. This is to turn on/off the red LED light and voice reminder of high temperature. Press **+** or **-** to turn on or turn off. After setting, press **SET** again or wait 5 seconds for it to finish automatically.

3-1. Temperature Measurement Sensitivity Value Adjustment : Three-click **SET** to High Temperature Warning icon, first press and hold **-** (without letting go) and then press **+**, then the sensitivity value (factory setting is 130) will be displayed on the screen and then let go. The user can adjust the value according to their needs.

3-2. **Increase Sensitivity** : Some materials of clothing may absorb near-infrared rays, which may cause component test delay or distance measurement error. If this occurs, the sensitivity setting value can be increased.

Press **-** to decrease the value, **the lower the value, the higher the sensitivity. The higher the sensitivity, the longer the measurement distance is allowed**, and vice versa. After setting, press **SET** again or wait 5 seconds for it to finish automatically. During the adjustment process, user can carry out temperature measurement to see whether it conforms to the user environment until it is satisfied. ***It is recommend that the subject stand at 50cm from the machine to avoid temperature measurement errors.**

3-3. **Decrease Sensitivity** : Press **+** to increase the value, **the higher the value, the lower the sensitivity**. After setting, press **SET** again or wait 5 seconds for it to finish automatically.


4. Multi-lingual Voice Prompts : Four-click **SET** to Multi-lingual Voice Prompts. When L 1 appears on the screen, press **+** or **-** to select language. The default language order is:

English (L1), Chinese (L2), French (L3), German (L4), Japanese (L5), Spanish (L6), Portuguese (L7), and Korean (L8). The default recognition voice is 'Hello' . After setting, press **SET** again or wait 5 seconds for it to finish automatically.

5. Volume : Five-click **SET** to Volume icon. Press **+** or **-** to adjust the volume level, 1 is mute, 5 is the loudest. After setting, press **SET** again or wait 5 seconds for it to finish automatically.

6. DC Mode : Six-click **SET** to Power icon. This is the symbol that appears when the device uses DC power. **(Please prepare your own USB charger above 5V_1A.)**

7. Battery Mode : Seven-click **SET** to Battery icon. This is the symbol that appears when the device uses batteries. If the battery voltage <4.6V, the battery symbol will flash, and it is recommended to replace the battery. It will automatically shut down after a beep sound (1 second) when the battery voltage <4.2V.

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